


Parts of the Dissecting Microscope




- A - Eyepiece lens, 10x
- B - Focus knob
- C - ON/OFF switch
- D - Objective lens, 2X or 4X
- E - Light
- F - Stage

How to use a Dissecting Microscope

1


Use the On/Off switch to turn the light source on.



2

Adjust the base

Depending on the sample, it may be better to use a black or a white background. Try and see which one is better.



3

Place the sample directly under the microscope lens.



4

Adjust eyepieces

Move the eyepieces until they are a comfortable distance apart for you to look down.



5

Focus

Turn the focus knob until the sample looks clear.




6

Use your hands to move the sample until you see something you want to investigate.



7

You have two magnification settings, 2X and 4X. Turn the objective lens to increase or decrease magnification.



8

When you have finished using the microscope, turn off the light to prevent it from overheating.

